


DECONTAMINATION PROCESS

Cleaning and disinfection of equipment are essential to reduce the risk of infection.

- Single-patient use items may be used more than once by the same patient.
- Single use items  must be discarded after each use.
- Follow cleaning technique bellow:



Is the equipment contaminated with body fluids (i.e. blood, urine, faeces, vomitus or phlegm)?

Was the equipment used for bathing/showering or toileting the patient?

Yes

No

Can the equipment be laundered?

Yes

RED alginate bag

Tie at the point of care and place in blue Elis bag for



No

Use Tristle solution to disinfect equipment.



Please complete the "I am Clean" sticker. Signed and dated.



Clean the equipment with **Clinell wipes**



Allow to dry before using it again or before storing.